



Instructions for Making Surgery Caps for Greenville Hospital System

Use only 100% washed cotton.

Cut circles. Use a serger with a differential feed to gather the edges of the fabric, or sew a basting stitch 1/4 inch in from edge and gather slightly, allowing fabric to turn to make a casing to hold 1/4 inch elastic. Stitch casing leaving a small opening to thread elastic into casing.

Circle Sizes and Elastic Lengths

Small: 16 inch circle	18 to 24 inches of elastic
Medium: 18 inches	22 to 24 inches of elastic
Large: 21 inches	26 to 28 inches of elastic